

MEDesign® Patient Transfer Belt

Product Reference: PTB



Description

The MEDesign® Patient Transfer Belt is a webbing belt with a two part, quick release buckle. It is an inexpensive means of assisting the carer to help a patient to move without grasping their clothing. It allows for a more secure grip and brings the patient closer.

Use

Place the belt around the patient's waist and snap close the buckle - ensure buckle is fully engaged before use. Adjust to a snug fit (do not over tighten) by pulling on the free strap which will lock in place and remain secure at the desired length. To release the belt, squeeze on the sides of the buckle and it will disengage.

Specification

Webbing: Black polypropylene webbing
Length 138cm (54") x Width 5cm (2")
Buckle: 50mm black acetal two part side release
Weight: 85g (3oz)



Safety Checklist

Ensure that webbing is not frayed. Ensure belt is attached securely. Ensure buckle is undamaged. Ensure that equipment is suitable for the client. Remove from service immediately if there is any sign of damage or fraying of the webbing. If in doubt do not use.

Care & Cleaning

To clean the belt, wash in warm water with an approved detergent, rinse well and hang up to dry. Clean the belt on a regular basis or after known contamination. Do not soak in disinfecting media, do not Autoclave, do not irradiate the Patient Transfer Belt.

Caution

Not suitable for use during post operative period. The risks associated with moving and handling are complex. Each situation must be assessed on its own merits. It is unwise for carers just to follow the instructions in our literature without proper assessment of the individual case. MEDesign® cannot accept responsibility for any consequences which may result from decisions made solely upon the basis of advice given herein.

©, MEDesign Ltd, 2002

Issue C12

