

MEDesign® Deluxe Reversible Prism Glasses

Product Code: BGD



Description

Our deluxe prism glasses have reversible lenses so the user can choose whether to look up or down. Read or watch TV whilst lying down. Helps alleviate strain in the neck and back. Used for belaying while rock climbing. Have also been used in MRI¹ scanners to alleviate the effects of claustrophobia. Can be worn with most normal spectacles.

1 Always check suitability/compatibility before use

Use

Place the prism glasses on the head just like a standard pair of spectacles and look straight ahead. They can be worn over a normal pair of spectacles if required. The lenses can easily be reversed so the wearer can choose whether to look up or down - see adjustment instructions on rear page.



Specification

Prism: Manufactured to optical standards

Weight: 57g (2oz)



Care & Cleaning

Treat the reversible prism glasses like a normal pair of spectacles. Wash in warm water with a mild detergent, rinse well and dry thoroughly. The prism lenses can be cleaned with a standard spectacle lens cleaning cloth. Do not soak in disinfecting media, do not Autoclave, do not irradiate the prism glasses. Inspect the glasses before each use to ensure they are whole, undamaged and safe to use.

©, MEDesign Ltd, 1982-2016

Issue C26



Ask for a copy of the MEDesign Back Pack catalogue showing our full product range.
MEDesign Ltd, Clock Tower Works, Railway Street, Southport, Merseyside, PR8 5BB, England.

Tel: +44 (0)1704 542373 Fax: +44 (0)1704 545214 Website: www.medesign.co.uk E-Mail: sales@medesign.co.uk



Reversible Lens Adjustment Instructions


Please note: MEDesign's Deluxe Reversible Prism Glasses are not a toy and are intended to be used by responsible adults. Care must be taken when releasing and/or reinserting the two prism holders from their respective locating slots in the rim of the spectacle frame.

Always ensure that any lens being removed cannot drop or fall onto any surface that could result in damage. Working carefully and methodically remove one prism lens at a time from the two locating slots in the rim of the spectacle frame. Use gentle uniform pressure to carefully and slowly push/slide the lens and its holder directly upwards in the two locating slots.

All Reversible Prism Glasses are closely inspected and quality tested during manufacture and then again prior to despatch to ensure satisfactory release and insertion of the prism lens holders from the spectacle frame.

Never try and force a lens in or out of the locating slots since this may result in damage to the prism holder or spectacle frame which is not covered by guarantee. If during use, or as a result of cleaning, the lens holder appears to stick or bind in the locating slot then we advise applying a very small amount of a suitable spectacle lubricant to the two side locating lug arms on the prism lens holder.

Locating lug arm
lubrication point



Locating lug arm
lubrication point

©, MEDesign Ltd, 1982-2016

Issue C26

